

- distributors
- technical information
- free samples
- new developments



Alconox®

Powdered Precision Cleaner

- Concentrated to save you money
- Biodegradable and readily disposable
- Penetrating wetting power to save you time
- Replaces corrosive acids and hazardous solvents
- Free rinsing to give you reliable results and no interfering residues
- Use to pass your cleaning validation tests for lab accreditation and plant inspection approval

Used to clean: Healthcare instruments, laboratory ware, vacuum equipment, tissue culture ware, personal protective equipment, sampling apparatus, catheters, tubing, pipes, radioactive contaminated articles, optical parts, electronic components, pharmaceutical apparatus, cosmetics manufacturing equipment, metal castings, forgings and stampings, industrial parts, tanks and reactors. Authorized by USDA for use in federally inspected meat and poultry plants. Passes inhibitory residue test for water analysis. FDA certified.

Used to remove: Soil, grit, grime, buffing compound, slime, grease, oils, blood, tissue, salts, deposits, particulates, solvents, chemicals, radioisotopes, radioactive contaminations, silicon oils, mold release agents.

Surfaces cleaned: Corrosion inhibited formulation recommended for glass, metal, stainless steel, porcelain, ceramic, plastic, rubber and fiberglass. Can be used on soft metals such as copper, aluminum, zinc and magnesium if rinsed promptly. Corrosion testing may be advisable.

Cleaning method: Soak, brush, sponge, cloth, ultrasonic, flow through clean-in-place. Will foam—not for spray or machine use.

Directions: Make a fresh 1% solution (2 1/2 Tbsp. per gal., 1 1/4 oz. per gal. or 10 grams per liter) in cold, warm, or hot water. If available use warm water. Use cold water for blood stains. For difficult soils, raise water temperature and use more detergent. Clean by soak, circulate, wipe, or ultrasonic method. Not for spray machines, will foam. For nonabrasive scouring, make paste. Use 2% solution to soak frozen stopcocks. To remove silver tarnish, soak in 1% solution in aluminum container. RINSE THOROUGHLY—preferably with running water. For critical cleaning, do final or all rinsing in distilled, deionized, or purified water. For food contact surfaces, rinse with potable water. Used on a wide range of glass, ceramic, plastic, and metal surfaces. Corrosion testing may be advisable.

Convenient Sizes:	Alconox Cat. #
Case 9 x 4 lb. Boxes	1104
25 lb. Carton	1125
50 lb. Carton	1150
100 lb. Drum	1101
300 lb. Drum	1103
Case 12 Box 50 x 1/2 oz.	1112
1 lb. makes 13 gal. cleaning solution	



Alconox is available from leading laboratory, hospital, clinical and industrial suppliers. To find a distributor for Alconox, Inc. detergents, visit "Find Dealer" at the website. To request FREE samples visit our Sample Request at www.alconox.com, write or call Alconox, Inc. today.

- distributors
- technical information
- free samples
- new developments

PHYSICAL DATA

	TYPICAL VALUE
pH of 1% solution	9.5
Flash Point (degrees F)	None
Phosphate Content (as Phosphorus)	7.3%
Organic Carbon (% calculated w/w)	11%
Fragrance Content	0%
Surface Tension 1% Sol'n (Dyne/cm)	32
Percent active ingredients	100%
Color: White and cream colored flakes	
Form: Powder	
Solubility in Water: To 10% (w/w) at ambient temperature	
Hard Water Effectiveness: Highly Effective	
Biodegradability: Biodegradable	
Foam Tendency: High Foaming	
Shelf Life: Two years from the date of manufacture	

Chemical Description

Alconox consists primarily of a homogeneous blend of sodium linear alkylaryl sulfonate, alcohol sulfate, phosphates, carbonates. Alconox is anionic in nature.

Cleaning Validation Methods:

Test a parameter of rinse water before and after rinsing the cleaned surface. No significant change in the parameter indicates no detectable residue. Parameters measured include: pH, conductivity, UV, TOC, HPLC, sodium concentration, phosphorus concentration, anionic surfactant concentration using inexpensive detergent water testing kits, surface tension, and surface analysis. Pharmaceutical Cleaning Validation references are available in the tech info section of www.alconox.com.

Health Safety Information:

OSHA Hazardous Ingredients: None

RCRA Hazard Class: Non-hazardous

Flammability: Non-flammable

Latex Content: None in detergent, packaging materials or adhesives.

Oral Toxicity: LD₅₀ > 5000 mg/kg oral rat

No ingredient defined as an oral toxicant by OSHA

Eye Irritation: Mild to Moderate eye irritant if not rinsed

Inhalation Toxicity: Non-irritating solution, powder a potential irritant

VOC Content: 0%

Carcinogenicity:

NTP = No IARC = No OSHA = No

All ingredients in Alconox are listed in TSCA inventory. USDA NSF cat A1

Precautions:

No special precautions other than good industrial hygiene and safety practices employed with any industrial chemical (see Directions). A Material Safety Data Sheet is available at www.alconox.com or by calling fax-on-demand at 914-948-4040 and following the prompts. Use fax document # 0311.

Contact Alconox, Inc. for purchase specifications. Typical data is not a specification.

While the information in this report should not be considered to be a product warranty, we urge you to investigate, test and verify the suitability of Alconox detergents for your specific application. We, of course, can not give permission to use, or recommend the use of, our detergents where they infringe patents. No representation or warranty is made as to the safety of products or materials mentioned under the Federal Food Additives Amendment of 1958.