

24-Hour Emergency Number  
(ChemTel) 888.255.3924 in U.S.

[cleaning@alconox.com](mailto:cleaning@alconox.com)

[alconox.com](http://alconox.com)

- distributors
- technical information
- free samples
- new developments
- white papers
- prices



## LUMINOX® Low-Foaming Neutral Cleaner

- Low-foaming for use in mechanical and spray washers
- Neutral pH for non-corrosive cleaning, safe handling and easy waste disposal
- No chelation or alkalinity waste treatment problems
- Concentrated for economical use
- Free rinsing to eliminate interfering residues
- Biodegradable and phosphate free

**Used to clean:** Aluminum and soft metal surfaces, medical instruments in aluminum trays, alkaline-sensitive lab and glassware, sensitive aerospace materials, pharmaceutical process equipment, medical devices, electronic components. FDA certified.

**Used to remove:** Oils, dirt, grime, adhesives and uncured epoxy, particulates, residues, fluxes, chemical solvents, machining oils.

**Surfaces cleaned:** Corrosion inhibited formulation recommended for glass, ferrous metals, non-ferrous metals, aluminum, copper, zinc, and other soft metals, ceramic, plastic, circuit boards and electronic components.

**Cleaning method:** Soak tanks, manual cleaning, ultrasonic tanks, mechanical washers, washer/sanitizers, spray washers, CIP systems.

**Directions:** For manual, soak and ultrasonic use, make a 3–5% solution (4–6 oz per gal, or 30–50 ml per L) in warm or hot water. Protective gloves and eyewear are recommended. For spray machines, add detergent at the hot water wash cycle to make a 2% solution (2 1/2 oz per gal, or 20 ml per L). For difficult soils raise water temperature

and use more detergent. RINSE THOROUGHLY—preferably with running water. For critical cleaning do final or all rinsing in distilled, deionized or purified water. Corrosion testing may be advisable. Material compatibility guidelines available on request.

Convenient Sizes	Alconox Cat. #
5 gal (19L) jerrycan	1905
Case of 4 x 1 gal (4 x 3.8L)	1901
15 gal (57L) drum	1915
55 gal (208L) drum	1955
1 gal makes 20–50 gal (1 L makes 20–50L) solution	



LUMINOX is available from leading laboratory, hospital, clinical and industrial suppliers. To find a distributor for Alconox, Inc. detergents, use the Locate Dealers search tool in the Support section of [alconox.com](http://alconox.com). To request FREE samples, use the Get Samples form at [alconox.com](http://alconox.com). Or call Alconox, Inc. at +914.948.4040, today.

LUTB.1

24-Hour Emergency Number  
(ChemTel) 888.255.3924 in U.S.  
cleaning@alconox.com

alconox.com

- distributors
- technical information
- free samples
- new developments
- white papers
- prices

Physical Data	Typical Value
<b>pH</b> (as is)	7.0
<b>Surface Tension 1% Sol'n</b> (dyne/cm)	34
<b>Conductivity 1% Sol'n</b> (mS/cm at 21 degrees C)	0.4
<b>Form</b>	Liquid
<b>Color</b>	Clear
<b>Solubility in Water</b>	Completely soluble in all proportions
<b>Hard Water Effectiveness</b>	Highly effective
<b>Biodegradability</b>	Yes, meets EC648/2004
<b>Foam Tendency</b>	Low foaming
<b>Shelf Life</b>	24 months from end of month of manufacture
<b>Fragrance and Dye Content</b>	0%
<b>Phosphate Content</b> (as phosphorus)	0%
<b>Flash Point</b> (degrees F)	>200 degrees F
<b>Density</b> (lb/gal)	8.75

### Chemical Description

LUMINOX is an aqueous, low-ionic content, phosphate-free blend of organic solvents, surfactants, and rinse aids.

### Cleaning Validation Methods

To quantify trace levels of Luminnox detergent for cleaning validation purposes, use HPLC. For details see the Cleaning Validation page in Support at alconox.com.

### Health Safety Information

<b>OSHA Hazardous Ingredients</b>	None
<b>RCRA Hazard Class</b>	Nonhazardous
<b>EPA Priority Pollutants</b>	None
<b>DOT Hazard Class</b>	Nonhazardous
<b>Flammability</b>	Nonflammable
<b>Latex Content</b>	None in product or packaging.
<b>Irritation</b>	Cat. 2 skin, Cat. 2A eye irritant
<b>Oral Toxicity</b>	LD <sub>50</sub> > 5000 mg/kg oral rat.
<b>Inhalation Toxicity</b>	Nonirritating
<b>VOC Content</b>	Less than 50 g/liter even at 10% use concentration.
<b>GMO Content</b>	None
<b>BSE/TSE Source</b>	None
<b>2007/68/EC Annex IIIa Allergens</b>	None
<b>REACH Art. 57 (SVHC)</b>	None
<b>RoHS</b>	No RoHS hazardous ingredients.
<b>Carcinogenicity</b>	<b>NTP</b> = No <b>IARC</b> = No <b>OSHA</b> = No
<b>Ingredients on TSCA Inventory</b>	Yes
<b>EC 907/2006</b>	Contains no enzymes, optical brighteners, perfumes, allergenic fragrances, or preservative agents.

### Precautions

No special precautions other than good industrial hygiene and safety practices employed with any industrial chemical (see Directions). A Safety Data Sheet is available at alconox.com.

Contact Alconox, Inc. for purchase specifications. Information presented is typical. Not to be taken as specifications. Typical data is not a specification.

While the information in this report should not be considered to be a product warranty, we urge you to investigate, test and verify the suitability of Alconox, Inc. detergents for your specific application. We, of course, can not give permission to use, or recommend the use of, our detergents where they infringe patents. No representation or warranty is made as to the safety of products or materials mentioned under the Federal Food Additives Amendment of 1958.

LUTB.1