

24-Hour Emergency Number  
(ChemTel) 888.255.3924 in U.S.

[cleaning@alconox.com](mailto:cleaning@alconox.com)

[alconox.com](http://alconox.com)

- distributors
- technical information
- free samples
- new developments
- white papers
- prices



## LIQUINOX® Critical-Cleaning Liquid Detergent

- Concentrated to save you money
- Replaces corrosive acids and hazardous solvents
- Phosphate-free, biodegradable and readily disposable
- Free rinsing to provide reliable results and no interfering residues
- For FDA cleaning validations and food contact surface cleaning

**Used to clean:** Healthcare instruments, laboratory ware, vacuum equipment, tissue culture ware, personal protective equipment, sampling apparatus, catheters, tubing, wine glasses, clean rooms, medical devices, optical parts, electronic components, pharmaceutical and cosmetics manufacturing equipment, metal castings, forgings and stampings, industrial parts, pipes, tanks and reactors. Used for phosphate sensitive analysis ware.

**Used to remove:** Soils, grit, grime, slime, grease, oils, blood, tissue, particulates, deposits, chemicals and solvents.

**Surfaces cleaned:** Corrosion inhibited formulation recommended for glass, metal, stainless steel, porcelain, ceramic, plastic, cement and fiber glass. Can be used on soft metals such as copper, aluminum, zinc and magnesium if rinsed promptly. Used for art restoration. Corrosion testing may be advisable.

**Cleaning method:** Soak, brush, sponge, cloth, ultrasonic, immersion/ recirculation cleaning in place. Will foam—not for spray or machine use.

**Directions:** Make a fresh 1% solution (2 1/2 tsp per gal, 1 1/4 oz per gal, or 10 ml per L) in cold, warm or hot water. If available, use warm water. Use cold water for blood stains. For difficult soils, raise water temperature and use more detergent. Clean by soak, circulate,

wipe or ultrasonic method. Not for spray machines: will foam. RINSE THOROUGHLY—preferably with running water. For critical cleaning, do final or all rinsing in distilled, deionized or purified water. For food contact surfaces, rinse with potable water. Used on a wide range of glass, ceramic, plastic and metal surfaces. Corrosion testing may be advisable. Material compatibility guidelines available on request.

Convenient Sizes	Alconox Cat. #
Case 12 x 1 qt (12 x 0.95L)	1232
Case of 4 x 1 gal (4 x 3.8L)	1201
5 gal (19L) jerrycan	1205
15 gal (57L) drum	1215
55 gal (208L) drum	1255
1 gal of concentrate makes 100 gal (1 L makes 100L) solution	



LIQUINOX passes inhibitory residue test for water analysis and is FDA certified. LIQUINOX is available from leading laboratory, hospital, clinical and industrial suppliers. To find a distributor for Alconox Inc. detergents, use the Locate Dealers search tool in the Support section of [alconox.com](http://alconox.com). To request FREE samples, use the Get Samples form at [alconox.com](http://alconox.com). Or call Alconox Inc. at +914.948.4040, today.

LQTB.3

24-Hour Emergency Number  
(ChemTel) 888.255.3924 in U.S.

[cleaning@alconox.com](mailto:cleaning@alconox.com)

[alconox.com](http://alconox.com)

- distributors
- technical information
- free samples
- new developments
- white papers
- prices

LIQUINOX Typical Data	Typical Value
<b>pH</b> (as is)	8.5
<b>Surface Tension 1% Sol'n</b> (dyn/cm)	29
<b>Conductivity 1% Sol'n</b> (mS/cm at 21 degrees C)	0.7
<b>Form</b>	Liquid
<b>Color</b>	Pale yellow
<b>Solubility in Water</b>	Completely soluble in all proportions
<b>Hard Water Effectiveness</b>	Highly effective
<b>Biodegradability</b>	Yes, meets EC648/2004
<b>Foam Tendency</b>	High foaming
<b>Shelf Life</b>	24 months from end of month of manufacture
<b>Fragrance and Dye Content</b>	0%
<b>Phosphate Content</b> (as phosphorus)	0%
<b>Flash Point</b> (degrees F)	None
<b>Organic Carbon</b> (% calculated w/w)	19%
<b>Density</b> (lb/gal)	8.9

### Chemical Description

Liquinox consists primarily of a homogeneous blend of sodium linear alkylaryl sulfonate, sodium xylene sulfonate, and ethoxylated alcohol. Liquinox contains a mixture of nonionic and anionic surfactants.

### Cleaning Validation and Residue Detection

To quantify trace levels of Liquinox for cleaning validation purposes, use conductivity, UV, TOC, or HPLC. For details see the Cleaning Validation page in Support at [alconox.com](http://alconox.com). For lab accreditation tests to show the detergent is rinsed away, see the Accreditation Residue Detection in the Support section of [alconox.com](http://alconox.com).

### Health and Safety Information

<b>OSHA Hazardous Ingredients</b>	None
<b>RCRA Hazard Class</b>	Nonhazardous
<b>EPA Priority Pollutants</b>	None
<b>DOT Hazard Class</b>	Nonhazardous
<b>Flammability</b>	Nonflammable
<b>Latex Content</b>	None in product or packaging.
<b>Irritation</b>	Cat. 2 skin, Cat. 2A eye irritant
<b>Oral Toxicity</b>	LD <sub>50</sub> > 5000 mg/kg oral rat.
<b>Inhalation Toxicity</b>	Nonirritating
<b>VOC Content</b>	0%
<b>GMO Content</b>	None
<b>BSE/TSE Source</b>	None
<b>2007/68/EC Annex IIIa Allergens</b>	None
<b>REACH Art. 57 (SVHC)</b>	None
<b>RoHS</b>	No RoHS hazardous ingredients.
<b>Carcinogenicity</b>	<b>NTP = No IARC = No OSHA = No</b>
<b>Ingredients on TSCA Inventory</b>	Yes
<b>EC 907/2006</b>	Contains no enzymes, optical brighteners, perfumes, allergenic fragrances, or preservative agents.

### Precautions

Avoid contact with skin, eyes and clothing. Do not take internally. In case of contact immediately flush skin with plenty of water. For eyes, flush with plenty of water for at least 15 minutes and get medical attention immediately. A Safety Data Sheet is available at [alconox.com](http://alconox.com). Contact Alconox Inc. for purchase specifications. Information presented is typical. Not to be taken as specifications. Typical data is not a specification.

While the information in this report should not be considered to be a product warranty, we urge you to investigate, test and verify the suitability of Alconox Inc. detergents for your specific application. We, of course, can not give permission to use, or recommend the use of, our detergents where they infringe patents. No representation or warranty is made as to the safety of products or materials mentioned under the Federal Food Additives Amendment of 1958.

LQTB.3