

# Alconox, Inc.

---

Critical cleaning detergents for laboratory, healthcare and industrial applications

30 Glenn Street, Suite 309 White Plains, NY 10603 USA

phone (914) 948-4040 fax (914) 948-4088 e-mail [cleaning@alconox.com](mailto:cleaning@alconox.com) website [www.alconox.com](http://www.alconox.com)

## ALCOJET® SAMPLE TRACE ANALYSIS UNDILUTED CONCENTRATE

<u>ELEMENT</u>	<u>DETECTED</u>	<u>DETECTION LIMIT</u>
Chlorine	7079 PPM	LESS THAN 100 PPM
Fluorine	NONE	LESS THAN 10 PPM
Bromine	NONE	LESS THAN 10 PPM
Iodine	NONE	LESS THAN 10 PPM
Boron	NONE	LESS THAN 10 PPM
Lead	NONE	LESS THAN 10 PPM
Zinc	NONE	LESS THAN 10 PPM
Cadmium	NONE	LESS THAN 50 PPM
Arsenic	NONE	LESS THAN 50 PPM
Antimony	NONE	LESS THAN 50 PPM
Bismuth	NONE	LESS THAN 10 PPM
Tin	NONE	LESS THAN 10 PPM
Sulfur	41 PPM	LESS THAN 10 PPM
Mercury	NONE	LESS THAN 0.5 PPM
Sodium	an ingredient in ALCOJET present at 20-30% by weight	
Phosphorus	an ingredient in ALCOJET present at 8.7% by weight	

Note: ALCOJET detergent is recommended for use at 1% dilutions.

