



30 Glenn Street, White Plains, NY 10603 USA

phone (914) 948-4040 fax (914) 948-4088 e-mail cleaning@alconox.com website www.alconox.com

DETERGENT 8®

Low Foaming Phosphate Free Detergent

SAMPLE TRACE ANALYSIS UNDILUTED CONCENTRATE

<u>ELEMENT</u>	<u>DETECTED</u>	<u>DETECTION LIMIT</u>
CHLORIDE	0.05 PPM	LESS THAN 0.01 PPM
FLUORIDE	None	LESS THAN 0.5 PPM
BROMIDE	0.09 PPM	LESS THAN 0.01 PPM
IODIDE	None	LESS THAN 1.0 PPM
SULFUR (as Sulfate)	None	LESS THAN 1.0 PPM
CHROMIUM	None	LESS THAN 1.0 PPM
HEXAVALENT CHROMIUM	None	LESS THAN 0.1 PPM
MERCURY	None	LESS THAN 0.1 PPM
LEAD	None	LESS THAN 0.05 PPM
ANTIMONY	None	LESS THAN 2.0 PPM
COPPER	None	LESS THAN 2.0 PPM
ZINC	None	LESS THAN 15.0 PPM
BISMUTH	None	LESS THAN 0.01 PPM
TIN	None	LESS THAN 1.0 PPM
CADMIUM	None	LESS THAN 0.25 PPM
ALUMINUM	None	LESS THAN 3.0 PPM
PHOSPHORUS	None	LESS THAN 20.0 PPM
SODIUM	7 PPM	LESS THAN 1.0 PPM
POTASSIUM	None	LESS THAN 1.0 PPM
SILICATES	None	LESS THAN 10 PPM
GALLIUM	None	LESS THAN 0.5 PPM
INDIUM	None	LESS THAN 0.5 PPM
SILVER	None	LESS THAN 1.0 PPM
NITRITE	None	LESS THAN 0.5 PPM
NITRATE	None	LESS THAN 5 PPM

Recommended dilution is 2-5% (3.5 oz. **Detergent 8** per gallon of water). When diluted, the amount of the above chemical traces will be diluted.

Detergent 8 is available in 1, 15, 30 and 55 gallons sizes
Catalog numbers 1701, 1715, 1730 and 1755

Malcolm McLaughlin,
VP Product and Business Development

