

Alconox, Inc.

Critical cleaning detergents for laboratory, healthcare and industrial applications

30 Glenn Street, Suite 309 White Plains, NY 10603 USA

phone (914) 948-4040 fax (914) 948-4088 e-mail cleaning@alconox.com website www.alconox.com

DET-O-JET® SAMPLE TRACE ANALYSIS UNDILUTED CONCENTRATE

<u>ELEMENT</u>	<u>DETECTED</u>	<u>DETECTION LIMIT</u>
Chlorine	9095 PPM	LESS THAN 100 PPM
Fluorine	NONE	LESS THAN 10 PPM
Bromine	NONE	LESS THAN 10 PPM
Iodine	NONE	LESS THAN 10 PPM
Boron	NONE	LESS THAN 10 PPM
Lead	NONE	LESS THAN 10 PPM
Zinc	NONE	LESS THAN 10 PPM
Cadmium	NONE	LESS THAN 50 PPM
Arsenic	NONE	LESS THAN 50 PPM
Antimony	NONE	LESS THAN 50 PPM
Bismuth	NONE	LESS THAN 10 PPM
Tin	NONE	LESS THAN 10 PPM
Sulfur	71 PPM	LESS THAN 10 PPM
Mercury	NONE	LESS THAN 0.5 PPM
Potassium	an ingredient in DET-O-JET present at 1-5% by weight	
Sodium	an ingredient in DET-O-JET present at 1-5% by weight	
Phosphorus	an ingredient in DET-O-JET present at 2.3% by weight	

Note: DET-O-JET detergent is recommended for use at 1% dilutions.

