

# Alconox, Inc.

Critical cleaning detergents for laboratory, healthcare and industrial applications

30 Glenn Street, White Plains, NY 10603 USA

phone (914) 948-4040 fax (914) 948-4088 e-mail [cleaning@alconox.com](mailto:cleaning@alconox.com) website [www.alconox.com](http://www.alconox.com)

Document #621

## SOLUJET™ SAMPLE TRACE ANALYSIS UNDILUTED CONCENTRATE

| <u>ELEMENT</u> | <u>DETECTED</u>                                    | <u>DETECTION LIMIT</u> |
|----------------|--|------------------------|
| Chlorine       | NONE   | LESS THAN 100 PPM      |
| Fluorine       | NONE   | LESS THAN 10 PPM       |
| Bromine        | NONE   | LESS THAN 10 PPM       |
| Iodine         | NONE   | LESS THAN 10 PPM       |
| Boron          | NONE   | LESS THAN 10 PPM       |
| Lead           | NONE   | LESS THAN 10 PPM       |
| Zinc           | NONE   | LESS THAN 10 PPM       |
| Cadmium        | NONE   | LESS THAN 50 PPM       |
| Arsenic        | NONE   | LESS THAN 50 PPM       |
| Antimony       | NONE   | LESS THAN 50 PPM       |
| Bismuth        | NONE   | LESS THAN 10 PPM       |
| Tin            | NONE   | LESS THAN 10 PPM       |
| Sulfur         | 60 PPM   | LESS THAN 10 PPM       |
| Mercury        | NONE   | LESS THAN 0.5 PPM      |
| Phosphorous    | NONE   | LESS THAN 100 PPM      |
| Potassium      | an ingredient in SOLUJET present at 1-5% by weight |                        |
| Sodium         | an ingredient in SOLUJET present at 1-5% by weight |                        |

Note: SOLUJET detergent is recommended for use at 1% dilutions.

