

Alconox, Inc.

Critical cleaning detergents for laboratory, healthcare and industrial applications

30 Glenn Street, White Plains, NY 10603 USA

phone (914) 948-4040 fax (914) 948-4088 e-mail cleaning@alconox.com website www.alconox.com

TERGAJET® SAMPLE TRACE ANALYSIS UNDILUTED CONCENTRATE

<u>ELEMENT</u>	<u>DETECTED</u>	<u>DETECTION LIMIT</u>
Chlorine	NONE	LESS THAN 20 PPM
Fluorine	NONE	LESS THAN 10 PPM
Bromine	NONE	LESS THAN 10 PPM
Iodine	NONE	LESS THAN 10 PPM
Boron	NONE	LESS THAN 10 PPM
Lead	NONE	LESS THAN 10 PPM
Zinc	NONE	LESS THAN 10 PPM
Cadmium	NONE	LESS THAN 10 PPM
Arsenic	NONE	LESS THAN 10 PPM
Antimony	NONE	LESS THAN 10 PPM
Bismuth	NONE	LESS THAN 10 PPM
Tin	NONE	LESS THAN 10 PPM
Mercury	NONE	LESS THAN 0.5 PPM
Phosphorus	NONE	LESS THAN 50 PPM
Sodium	an ingredient in TERGAJET present at 20-30% by weight	
Sulfur	sulfates are an ingredient in Tergajet at 5-15% by weight	

Note: TERGAJET detergent is recommended for use at 1% dilutions.

