



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 - Fax.914.948.4088

www.alconox.com - cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 10/3/2019

CERTIFICATE OF ANALYSIS
ALCONOX

Manufactured by: Alconox, Inc.

Sub lot No: A1K9 A1K9A, A1K9B, A1K9C, A1K9D, A1K9E

Sub lot numbers of Master lot No: A2H8

Manuf. Cat. no.: 1101, 1103, 1104, 1104-1, 1112, 1112-1, 1125, 1150

Date Of Manufacture: 10/1/2019

Expires end of: 10/21

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

| <u>ANALYSIS</u> | <u>RESULT</u> | <u>ACCEPTED RANGE</u> |
|---|----------------------|--|
| Appearance: | OK | White powder with cream colored flakes |
| Solubility: (1% solution at 25 deg C, with agitation) | COMP | Complete |
| pH:(1% solution in DI water) | 9.2 | 9.0 - 10.0 |
| Surface Tension:(1% solution Dyne/cm) | 31.2 | 29.0 - 33.0 |

Laura Falla, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.