



Alconox, Inc.

The Leader in Critical-Cleaning
30 Glenn St., Suite 309, White Plains, NY 10603 USA
Tel.914.948.4040 • Fax.914.948.4088
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 1/11/2016

CERTIFICATE OF ANALYSIS **ALCONOX**

Manufactured by: Alconox, Inc.

Sub lot No: A2A6 A2A6A, A2A6B, A2A6C, A2A6D, A2A6E

Sub lot numbers of Master lot No: A4A5

Manuf. Cat. no.: 1101, 1103, 1104, 1104-1, 1112, 1112-1, 1125, 1150

Date Of Manufacture: 1/7/2016

Shelf Life Expiration Date (730 days): 1/6/2018

Storage Conditions: Closed upright and in a cool dry place (15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:	OK	White powder with cream colored flakes
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH:(1% solution in DI water)	9.4	9.0 - 10.0
Surface Tension:(1% solution Dyne/cm)	31.3	29.0 - 33.0

Malcolm McLaughlin, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.