



Alconox, Inc.

The Leader in Critical-Cleaning  
30 Glenn St., Suite 309, White Plains, NY 10603 USA  
Tel.914.948.4040 • Fax.914.948.4088  
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 4/24/2017

**CERTIFICATE OF ANALYSIS**  
**ALCONOX**

**Manufactured by: Alconox, Inc.**

Sub lot No: A3D7 A3D7A, A3D7B, A3D7C, A3D7D, A3D7E

Sub lot numbers of Master lot No: A1A7

Manuf. Cat. no.: 1101, 1103, 1104, 1104-1, 1112, 1112-1, 1125, 1150

Date Of Manufacture: 4/19/2017

Expires end of: 04/19

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

**ANALYSIS**

**RESULT**

**ACCEPTED RANGE**

Appearance:	OK	White powder with cream colored flakes
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH:(1% solution in DI water)	9.4	9.0 - 10.0
Surface Tension:(1% solution Dyne/cm)	30.8	29.0 - 33.0

Malcolm McLaughlin, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*