



Alconox, Inc.

The Leader in Critical-Cleaning  
30 Glenn St., Suite 309, White Plains, NY 10603 USA  
Tel.914.948.4040 • Fax.914.948.4088  
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 5/17/2017

## **CERTIFICATE OF ANALYSIS**

### **ALCONOX**

#### **Manufactured by: Alconox, Inc.**

Sub lot No: A3E7 A3E7A, A3E7B, A3E7C, A3E7D, A3E7E

Sub lot numbers of Master lot No: A1A7

Manuf. Cat. no.: 1101, 1103, 1104, 1104-1, 1112, 1112-1, 1125, 1150

Date Of Manufacture: 5/16/2017

Expires end of: 05/19

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

#### **ANALYSIS**

#### **RESULT**

#### **ACCEPTED RANGE**

Appearance:	OK	White powder with cream colored flakes
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH:(1% solution in DI water)	9.2	9.0 - 10.0
Surface Tension:(1% solution Dyne/cm)	30.6	29.0 - 33.0

Malcolm McLaughlin, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*