



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 - Fax.914.948.4088

www.alconox.com - cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 1/27/2020

**CERTIFICATE OF ANALYSIS**  
**CITRANOX**

**Manufactured by: Alconox, Inc.**

Sub lot No: C1A0 C1A0F, C1A0J, C1A0K, C1A0L, C1A0M

Sub lot numbers of Master lot No: C1F8

Manuf. Cat. no.: 1801, 1801-1, 1805, 1815, 1830, 1855

Date Of Manufacture: 1/27/2020

Expires end of: 01/22

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<b><u>ANALYSIS</u></b>	<b><u>RESULT</u></b>	<b><u>ACCEPTED RANGE</u></b>
Appearance:	OK	Clear, olive colored liquid
pH: (1% solution in DI water)	2.8	2.0 - 3.0
Surface Tension:(1% solution Dyne/cm)	32.3	29.0 - 33.0
Specific Gravity:	1.11	1.09 - 1.14

**Laura Falla, QA**

**Alconox, Inc.**

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*