



Alconox, Inc.

The Leader in Critical-Cleaning

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 11/9/2018

## **CERTIFICATE OF ANALYSIS**

### **CITRANOX**

**Manufactured by: Alconox, Inc.**

Sub lot No: C1L8 C1L8F, C1L8J, C1L8K, C1L8L, C1L8M

Sub lot numbers of Master lot No: C2K7

Manuf. Cat. no.: 1801, 1801-1, 1805, 1815, 1830, 1855

Date Of Manufacture: 11/7/2018

Expires end of: 11/20

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<u>ANALYSIS</u>	<u>RESULT</u>	<u>ACCEPTED RANGE</u>
Appearance:	OK	Clear, olive colored liquid
pH: (1% solution in DI water)	2.8	2.0 - 3.0
Surface Tension:(1% solution Dyne/cm)	32.1	29.0 - 33.0
Specific Gravity:	1.115	1.09 - 1.14

Malcolm McLaughlin, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*

COA NUMBER: 7648