



CERTIFICATE OF ANALYSIS
CITRANOX
FOR BECKMAN COULTER (BCI)

Date Of Manufacture: 01/07/2016
Manufactured by: Alconox, Inc.
Manufacturer Lot number: C9A6
BCI Lot number: 534862
Manuf. Cat. No.: 1801B
BCI Part No: 81912
Shelf Life Expiration Date: 12/31/2017 To meet BC requirements
Storage Conditions: Closed upright and in a cool dry place (15 - 30 deg C or 60 - 90 deg F).

<u>Value</u>	<u>-</u>	<u>Result</u>	<u>Acceptance Range</u>
Appearance:		OK	Clear pale olive colored liquid
pH: (1% solution in DI water)		2.7	2.0-3.0
Surface Tension: (1% solution Dyne/cm)		31.3	29.0 – 33.0
Specific Gravity		1.112	1.090 – 1.140
Number of Bar Code Labels Used for Labels		720	Not Applicable
Number of Bar Code Labels Used for Records		3	Not Applicable
Number of Bar Code Labels Discarded		19	Not Applicable



Malcolm McLaughlin, QA



CITRANOX

REF 81912

For use in the Access Immunoassay Systems.

 2017-12-31

LOT 534862



Alconox, Inc.
The Leader in Critical Cleaning
30 Glenn St. White Plains NY 10603 USA
Tel 914.948.4040 • Fax 914-948-4088
www.alconox.com • cleaning@alconox.com

Distributed by: BECKMAN COULTER, INC.
Brea, CA 92821 U.S.A.

Made In U.S.A. Printed In U.S.A.
©2015 Beckman Coulter, Inc. Beckman Coulter, Inc., Brea, CA 92821

101550A