



**CERTIFICATE OF ANALYSIS  
CITRANOX  
FOR BECKMAN COULTER (BCI)**

**Date Of Manufacture: 02/20/2018**  
**Manufactured by: Alconox, Inc.**  
**Manufacturer Lot number: C1B8**  
**BCI Lot number: C9C8**  
**Manuf. Cat. No.: 1801B**  
**BCI Part No: 81912**  
**Shelf Life Expiration Date: 01/31/20 to meet BC requirements**  
**Storage Conditions: Closed upright and in a cool dry place (15 - 30 deg C or 60 - 90 deg F).**

<u>Value</u>	<u>Result</u>	<u>Acceptance Range</u>
Appearance:	OK	Clear pale olive colored liquid
pH: (1% solution in DI water)	2.8	2.0-3.0
Surface Tension: (1% solution Dyne/cm)	31.6	29.0 - 33.0
Specific Gravity	1.116	1.090 - 1.140
Number of Bar Code Labels Used for Labels	720	Not Applicable
Number of Bar Code Labels Used for Records	3	Not Applicable
Number of Bar Code Labels Destroyed	483	Not Applicable

Malcolm McLaughlin, QA

**REF** 81912

**2020-01-31**

**LOT** C9C8



••H628819121A••



••552013120C9C8A7••



(01)15099590223571  
(17)200131  
(10)C9C8

**Alconox, Inc.**  
**The Leader in Critical Cleaning**  
30 Glenn St White Plains NY 10603 USA  
Tel 914.948.4040 • Fax 914-948-4088  
[www.alconox.com](http://www.alconox.com) • [cleaning@alconox.com](mailto:cleaning@alconox.com)