



CERTIFICATE OF ANALYSIS
CITRANOX
FOR BECKMAN COULTER (BCI)

Date Of Manufacture: 08/11/2016
Manufactured by: Alconox, Inc.
Manufacturer Lot number: C9H6
BCI Lot number: 630001
Manuf. Cat. No.: 1801B
BCI Part No: 81912
Shelf Life Expiration Date: 07/2018 To meet BC requirements
Storage Conditions: Closed upright and in a cool dry place (15 - 30 deg C or 60 - 90 deg F).

<u>Value</u>	<u>-</u>	<u>Result</u>	<u>Acceptance Range</u>
Appearance:		OK	Clear pale olive colored liquid
pH: (1% solution in DI water)		2.9	2.0-3.0
Surface Tension: (1% solution Dyne/cm)		31.1	29.0 – 33.0
Specific Gravity		1.106	1.090 – 1.140
Number of Bar Code Labels Used for Labels		720	Not Applicable
Number of Bar Code Labels Used for Records		3	Not Applicable
Number of Bar Code Labels Destroyed		32	Not Applicable

Malcolm McLaughlin, QA

REF 81912

2018-07-31

LOT 630001



+H628819121A



+\$\$2073118630001AW



(01)15099590223571

(17)180731

(10)630001