



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 3/5/2020

CERTIFICATE OF ANALYSIS
DETOJET

Manufactured by: Alconox, Inc.

Sub lot No: D1B0 D1B0F, D1B0G, D1B0H, D1B0J, D1B0K, D1B0L, D1B0M

Sub lot numbers of Master lot No: D1J8

Manuf. Cat. no.: 1601, 1601-1, 1605, 1615, 1630, 1632, 1632-1, 1655

Date Of Manufacture: 2/28/2020

Expires end of: 02/22

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:

OK

Clear, water white liquid

pH: (1% solution in DI water)

12.1

11.5 - 13.0

Foam (1% in DI water)

NIL

NIL

Specific Gravity:

1.33

1.27 - 1.34

Laura Falla, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.