



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 - Fax.914.948.4088

www.alconox.com - cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 8/27/2019

CERTIFICATE OF ANALYSIS **DETOJET**

Manufactured by: Alconox, Inc.

Sub lot No: D1H9 D1H9F, D1H9G, D1H9H, D1H9J, D1H9K, D1H9L, D1H9M

Sub lot numbers of Master lot No: D1J8

Manuf. Cat. no.: 1601, 1601-1, 1605, 1615, 1630, 1632, 1632-1, 1655

Date Of Manufacture: 8/26/2019

Expires end of: 08/21

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<u>ANALYSIS</u>	<u>RESULT</u>	<u>ACCEPTED RANGE</u>
Appearance:	OK	Clear, water white liquid
pH: (1% solution in DI water)	12.2	11.5 - 13.0
Foam (1% in DI water)	NIL	NIL
Specific Gravity:	1.33	1.27 - 1.34

Laura Falla, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.