



Alconox, Inc.

The Leader in Critical-Cleaning  
30 Glenn St., Suite 309, White Plains, NY 10603 USA  
Tel.914.948.4040 • Fax.914.948.4088  
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 2/2/2016

## ***CERTIFICATE OF ANALYSIS*** **DETERGENT 8**

**Manufactured by: Alconox, Inc.**

Sub lot No: E1A6 E1A6G, E1A6J, E1A6K, E1A6L, E1A6M

Sub lot numbers of Master lot No: E1L4

Manuf. Cat. no.: 1701, 1701-1, 1705, 1715, 1755

Date Of Manufacture: 1/26/2016

Shelf Life Expiration Date (730 days): 1/25/2018

Storage Conditions: Closed upright and in a cool dry place (15 - 30 deg C or 60 - 90 deg F).

<b><u>ANALYSIS</u></b>	<b><u>RESULT</u></b>	<b><u>ACCEPTED RANGE</u></b>
Appearance:	OK	Clear, water white liquid
pH: (1% solution in DI water)	11.0	10.5 - 11.8
Surface Tension:(1% solution Dyne/cm)	31.8	30.0 - 34.0
Foam (1% in DI water)	NIL	NIL
Specific Gravity:	0.990	0.97 - 1.01

**Malcolm McLaughlin, QA**

**Alconox, Inc.**

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*