



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 - Fax.914.948.4088

www.alconox.com - cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 3/5/2019

CERTIFICATE OF ANALYSIS ***DETERGENT 8***

Manufactured by: Alconox, Inc.

Sub lot No: E1C9 E1C9G, E1C9J, E1C9K, E1C9L, E1C9M

Sub lot numbers of Master lot No: E1C7

Manuf. Cat. no.: 1701, 1701-1, 1705, 1715, 1755

Date Of Manufacture: 3/4/2019

Expires end of: 03/21

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<u>ANALYSIS</u>	<u>RESULT</u>	<u>ACCEPTED RANGE</u>
Appearance:	OK	Clear, water white liquid
pH: (1% solution in DI water)	11.0	10.5 - 11.8
Surface Tension:(1% solution Dyne/cm)	32.9	30.0 - 34.0
Foam (1% in DI water)	NIL	NIL
Specific Gravity:	0.99	0.97 - 1.01

Laura Falla, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.