



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 - Fax.914.948.4088

www.alconox.com - cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 8/20/2019

## ***CERTIFICATE OF ANALYSIS*** ***DETERGENT 8***

**Manufactured by: Alconox, Inc.**

Sub lot No: E1H9 E1H9G, E1H9J, E1H9K, E1H9L, E1H9M

Sub lot numbers of Master lot No: E1F8

Manuf. Cat. no.: 1701, 1701-1, 1705, 1715, 1755

Date Of Manufacture: 8/14/2019

Expires end of: 08/21

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<b><u>ANALYSIS</u></b>	<b><u>RESULT</u></b>	<b><u>ACCEPTED RANGE</u></b>
Appearance:	OK	Clear, water white liquid
pH: (1% solution in DI water)	11.1	10.5 - 11.8
Surface Tension:(1% solution Dyne/cm)	33.1	30.0 - 34.0
Foam (1% in DI water)	NIL	NIL
Specific Gravity:	0.99	0.97 - 1.01

**Laura Falla, QA**

**Alconox, Inc.**

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*