



Alconox, Inc.

The Leader in Critical-Cleaning
30 Glenn St., Suite 309, White Plains, NY 10603 USA
Tel.914.948.4040 • Fax.914.948.4088
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 5/24/2016

CERTIFICATE OF ANALYSIS **TERGAZYME**

Manufactured by: Alconox, Inc.

Sub lot No: G1E6 G1E6A, G1E6B, G1E6C, G1E6D

Sub lot numbers of Master lot No: G1A5

Manuf. Cat. no.: 1301, 1303, 1304, 1304-1, 1325, 1350

Date Of Manufacture: 5/19/2016

Expires end of: 05/18

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

| | | |
|---|------|---|
| Appearance: | OK | White powder with cream colored flakes and small brown granules |
| Solubility: (1% solution at 25 deg C, with agitation) | COMP | Complete |
| pH:(1% solution in DI water) | 9.2 | 9.0 - 10.0 |
| Surface Tension:(1% solution Dyne/cm) | 30.2 | 29.0 - 33.0 |

Malcolm McLaughlin, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.