



CERTIFICATE OF ANALYSIS
HYDROPHARM

Date Of Manufacture: 01/29/2016

Manufactured by: Alconox, Inc.

Lot number: H032415

Manuf. Cat. no.: HYD-55, HYD-30

Shelf Life Expiration Date (1095 days): 03/24/2018

Storage Conditions: Closed upright and in a cool dry place (15 - 30 deg C or 60 - 90 deg F).

<u>Value</u>	-	<u>Result</u>	<u>Acceptance Range</u>
Appearance:		OK	Clear pale olive colored liquid
pH: (1% solution in DI water)		2.7	2.0-3.0
Surface Tension: (1% solution Dyne/cm)		35.3	31.5-36.0
Specific Gravity		1.056	1.01 – 1.15

Malcolm McLaughlin