



Alconox, Inc.

The Leader in Critical-Cleaning
30 Glenn St., Suite 309, White Plains, NY 10603 USA
Tel.914.948.4040 • Fax.914.948.4088
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 3/2/2016

CERTIFICATE OF ANALYSIS
ALCOJET

Manufactured by: Alconox, Inc.

Sub lot No: J1B6 J1B6A, J1B6B, J1B6C, J1B6D

Sub lot numbers of Master lot No: J1A5

Manuf. Cat. no.: 1401, 1403, 1404, 1404-1, 1425, 1450

Date Of Manufacture: 2/29/2016

Expires end of: 02/18

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

| | | |
|---|------|--------------|
| Appearance: | OK | White powder |
| Solubility: (1% solution at 25 deg C, with agitation) | COMP | Complete |
| pH:(1% solution in DI water) | 11.4 | 11.0 - 12.5 |
| Surface Tension:(1% solution Dyne/cm) | 32.2 | 31.0 - 36.0 |
| Foam (1% in DI water) | NIL | NIL |

Malcolm McLaughlin, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.