



Alconox, Inc.

The Leader in Critical-Cleaning

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 7/31/2018

CERTIFICATE OF ANALYSIS LIQUINOX

Manufactured by: Alconox, Inc.

Sub lot No: L2G8 L2G8F, L2G8H, L2G8J, L2G8K, L2G8L, L2G8M

Sub lot numbers of Master lot No: L1H7

Manuf. Cat. no.: 1201, 1201-1, 1205, 1215, 1230, 1232, 1232-1, 1255

Date Of Manufacture: 7/30/2018

Expires end of: 07/20

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:	OK	Clear golden liquid
pH: (as is)	8.4	8.0 - 9.5
Surface Tension:(1% solution Dyne/cm)	28.4	27.0 - 31.0
Specific Gravity:	1.074	1.06 - 1.09

Malcolm McLaughlin, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.