



Alconox, Inc.

The Leader in Critical-Cleaning  
30 Glenn St., Suite 309, White Plains, NY 10603 USA  
Tel.914.948.4040 • Fax.914.948.4088  
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 9/27/2016

## **CERTIFICATE OF ANALYSIS** **LIQUINOX**

**Manufactured by: Alconox, Inc.**

Sub lot No: L2J6 L2J6F, L2J6H, L2J6J, L2J6K, L2J6L, L2J6M

Sub lot numbers of Master lot No: L2M5

Manuf. Cat. no.: 1201, 1201-1, 1205, 1215, 1230, 1232, 1232-1, 1255

Date Of Manufacture: 9/23/2016

Expires end of: 09/18

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<b><u>ANALYSIS</u></b>	<b><u>RESULT</u></b>	<b><u>ACCEPTED RANGE</u></b>
Appearance:	OK	Clear golden liquid
pH: (as is)	8.4	8.0 - 9.5
Surface Tension:(1% solution Dyne/cm)	29.0	27.0 - 31.0
Specific Gravity:	1.070	1.06 - 1.09

Malcolm McLaughlin, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*