



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 - Fax.914.948.4088

www.alconox.com - cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 12/31/2019

## ***CERTIFICATE OF ANALYSIS*** ***LIQUINOX***

**Manufactured by: Alconox, Inc.**

Sub lot No: L2M9 L2M9F, L2M9H, L2M9J, L2M9K, L2M9L, L2M9M

Sub lot numbers of Master lot No: L1J8

Manuf. Cat. no.: 1201, 1201-1, 1205, 1215, 1230, 1232, 1232-1, 1255

Date Of Manufacture: 12/30/2019

Expires end of: 12/21

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<b><u>ANALYSIS</u></b>	<b><u>RESULT</u></b>	<b><u>ACCEPTED RANGE</u></b>
Appearance:	OK	Clear golden liquid
pH: (as is)	8.5	8.0 - 9.5
Surface Tension:(1% solution Dyne/cm)	31.0	27.0 - 31.0
Specific Gravity:	1.08	1.06 - 1.09

**Laura Falla, QA**

**Alconox, Inc.**

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*