



Alconox, Inc.

The Leader in Critical-Cleaning

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 1/25/2018

## CERTIFICATE OF ANALYSIS LIQUINOX

**Manufactured by: Alconox, Inc.**

Sub lot No: L3A8 L3A8F, L3A8H, L3A8J, L3A8K, L3A8L, L3A8M

Sub lot numbers of Master lot No: L3M6

Manuf. Cat. no.: 1201, 1201-1, 1205, 1215, 1230, 1232, 1232-1, 1255

Date Of Manufacture: 1/24/2018

Expires end of: 01/20

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

**ANALYSIS**

**RESULT**

**ACCEPTED RANGE**

Appearance:	OK	Clear golden liquid
pH: (as is)	8.5	8.0 - 9.5
Surface Tension:(1% solution Dyne/cm)	28.9	27.0 - 31.0
Specific Gravity:	1.076	1.06 - 1.09

Malcolm McLaughlin, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*