



Alconox, Inc.

The Leader in Critical-Cleaning
30 Glenn St., Suite 309, White Plains, NY 10603 USA
Tel.914.948.4040 • Fax.914.948.4088
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 2/24/2016

CERTIFICATE OF ANALYSIS

LIQUINOX

Manufactured by: Alconox, Inc.

Sub lot No: L3B6 L3B6F, L3B6H, L3B6J, L3B6K, L3B6L, L3B6M

Sub lot numbers of Master lot No: L2A5

Manuf. Cat. no.: 1201, 1201-1, 1205, 1215, 1230, 1232, 1232-1, 1255

Date Of Manufacture: 2/22/2016

Expires end of: 02/18

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

| | | |
|---------------------------------------|-------|---------------------|
| Appearance: | OK | Clear golden liquid |
| pH: (as is) | 8.4 | 8.0 - 9.5 |
| Surface Tension:(1% solution Dyne/cm) | 31.3 | 29.0 - 33.0 |
| Specific Gravity: | 1.080 | 1.06 - 1.09 |

Malcolm McLaughlin, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.