



Alconox, Inc.

The Leader in Critical-Cleaning
30 Glenn St., Suite 309, White Plains, NY 10603 USA
Tel.914.948.4040 • Fax.914.948.4088
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 2/4/2016

CERTIFICATE OF ANALYSIS

LUMINOX

Manufactured by: Alconox, Inc.

Sub lot No: M1A6 M1A6G, M1A6J, M1A6L, M1A6M

Sub lot numbers of Master lot No: M1L4

Manuf. Cat. no.: 1901, 1901-1, 1915, 1955

Date Of Manufacture: 1/28/2016

Shelf Life Expiration Date (730 days): 1/27/2018

Storage Conditions: Closed upright and in a cool dry place (15 - 30 deg C or 60 - 90 deg F).

<u>ANALYSIS</u>	<u>RESULT</u>	<u>ACCEPTED RANGE</u>
Appearance:	OK	Clear liquid
pH: (as is)	7.0	6.0 - 8.0
Surface Tension:(1% solution Dyne/cm)	34.0	32.0 - 38.0
Specific Gravity:	1.030	1.00 - 1.06
Foam (1% in DI water)	NIL	NIL

Malcolm McLaughlin, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.