

2/12/2016

CERTIFICATE OF ANALYSIS LQ9

Manufactured by: Alconox, Inc.

Lot No: Q1B6

Manuf. Cat. no.: LQ941, LQ905, LQ915, LQ930, and LQ955

Date of Manufacture: 2/09/2016

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.

Shelf Life Expiration Date: 02/2018

Storage Conditions: Closed upright and in a cool dry place (15 - 30 deg C or 60 - 90 deg F).

<u>Analysis</u>	<u>Result</u>	<u>Acceptance Range</u>
Appearance:	OK	Clear golden liquid
pH: (as is)	8.3	8.0-9.5
Surface Tension: (1% solution dyne/cm)	31.3	29.0 – 33.0
Specific Gravity:	1.080	1.060 – 1.090

Malcolm McLaughlin



Alconox, Inc