



Alconox, Inc.

The Leader in Critical-Cleaning
30 Glenn St., Suite 309, White Plains, NY 10603 USA
Tel.914.948.4040 • Fax.914.948.4088
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 4/15/2016

CERTIFICATE OF ANALYSIS
LQ9

Manufactured by: Alconox, Inc.

Sub lot No: Q1D6

Sub lot numbers of Master lot No: None

Manuf. Cat. no.: LQ941, LQ905, LQ915, LQ955

Date Of Manufacture: 4/4/2016

Expires end of: 04/18

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<u>ANALYSIS</u>	<u>RESULT</u>	<u>ACCEPTED RANGE</u>
Appearance:	OK	Clear golden liquid
pH: (as is)	8.5	8.0 - 9.5
Surface Tension:(1% solution Dyne/cm)	31.2	29.0 - 33.0
Specific Gravity:	1.080	1.06 - 1.09

Malcolm McLaughlin, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.