



Alconox, Inc.

The Leader in Critical-Cleaning
30 Glenn St., Suite 309, White Plains, NY 10603 USA
Tel.914.948.4040 • Fax.914.948.4088
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 5/16/2016

CERTIFICATE OF ANALYSIS

SOLUJET

Manufactured by: Alconox, Inc.

Sub lot No: S1E6 S1E6F, S1E6G, S1E6J, S1E6K, S1E6L, S1E6M

Sub lot numbers of Master lot No: S1A5

Manuf. Cat. no.: 2101, 2101-1, 2105, 2115, 2130, 2155

Date Of Manufacture: 5/11/2016

Expires end of: 05/18

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<u>ANALYSIS</u>	<u>RESULT</u>	<u>ACCEPTED RANGE</u>
Appearance:	OK	pale straw colored liquid
pH: (1% solution in DI water)	12.1	11.5 - 13.0
Surface Tension:(1% solution Dyne/cm)	32.0	31.5 - 36.5
Specific Gravity:	1.278	1.26 - 1.34

Malcolm McLaughlin, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.