



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 - Fax.914.948.4088

www.alconox.com - cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 11/12/2019

## ***CERTIFICATE OF ANALYSIS*** **TERGAJET**

**Manufactured by: Alconox, Inc.**

Sub lot No: T1L9 T1L9A, T1L9B, T1L9C, T1L9D

Sub lot numbers of Master lot No: T1G8

Manuf. Cat. no.: 2201, 2203, 2204, 2204-1, 2225, 2250

Date Of Manufacture: 11/9/2019

Expires end of: 11/21

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<b><u>ANALYSIS</u></b>	<b><u>RESULT</u></b>	<b><u>ACCEPTED RANGE</u></b>
Appearance:	OK	White powder
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH:(1% solution in DI water)	11.6	11.0 - 12.5
Surface Tension:(1% solution Dyne/cm)	32.7	31.0 – 36.0
Foam (1% in DI water)	NIL	NIL

**Laura Falla, QA**

**Alconox, Inc.**

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*