



VWR.com
1.800.932.5000

Distributed by: VWR International
Manufacturer: Alconox, Inc.
30 Glenn St.
White Plains, NY 10603 USA
Ph 914-948-4040
cleaning@alconox.com

November 15, 2016

CERTIFICATE OF ANALYSIS
LABTONE

Manufactured by: Alconox, Inc.
Lot No: V1H6
VWR Cat. No.: 21850-027, 89030-112 (4X8 lb bottle size) and 21850-080, 89030-114 (25 lb box size)
Date Of Manufacture: 08/18/2016
Shelf Life Expiration Date (2 years from manufacture): 08/18
Storage Conditions: Closed upright and in a cool dry place (15 – 30 deg C or 60 – 90 deg F).

<u>ANALYSIS</u>	<u>RESULT</u>	<u>ACCEPTED RANGE</u>
Appearance:	OK	White powder
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH: (1% solution in DI water)	9.8	9.5 – 11.5
Surface Tension (dynes/cm): (1% solution in DI water)	29.8	28.0 – 33.0

Malcolm McLaughlin