



Alconox, Inc.

The Leader in Critical-Cleaning
30 Glenn St., Suite 309, White Plains, NY 10603 USA
Tel.914.948.4040 • Fax.914.948.4088
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 6/7/2016

CERTIFICATE OF ANALYSIS

DETONOX

Manufactured by: Alconox, Inc.

Sub lot No: X060616

Sub lot numbers of Master lot No: None

Manuf. Cat. no.: 2301, 2305, 2315, 2355

Date Of Manufacture: 6/6/2016

Expires end of: 06/18

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:	OK	Clear pale yellow liquid –nearly odorless
pH: (1% solution in DI water)	10.1	9.5 - 11.5
Surface Tension:(1% solution Dyne/cm)	30.9	29.0 - 33.0
Specific Gravity:	1.202	1.155 -1.225

Malcolm McLaughlin, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.