



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 9/30/2020

CERTIFICATE OF ANALYSIS LIQUINOX

Manufactured by: Alconox, Inc.

Sub lot No: L2J0 L2J0F, L2J0H, L2J0J, L2J0K, L2J0L, L2J0M

Sub lot numbers of Master lot No: L1H9

Manuf. Cat. no.: 1201, 1201-1, 1205, 1215, 1230, 1232, 1232-1, 1255

Date Of Manufacture: 9/25/2020

Expires end of: 09/22

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<u>ANALYSIS</u>	<u>RESULT</u>	<u>ACCEPTED RANGE</u>
Appearance:	OK	Clear golden liquid
pH: (as is)	8.5	8.0 - 9.5
Surface Tension:(1% solution Dyne/cm)	29.4	27.0 - 31.0
Specific Gravity:	1.08	1.06 - 1.09

Laura Falla, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.