



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 7/21/2020

CERTIFICATE OF ANALYSIS SOLUJET

Manufactured by: Alconox, Inc.

Sub lot No: S2G0 S2G0F, S2G0G, S2G0J, S2G0K, S2G0L, S2G0M

Sub lot numbers of Master lot No: S1H9

Manuf. Cat. no.: 2101, 2101-1, 2105, 2115, 2130, 2155

Date Of Manufacture: 7/20/2020

Expires end of: 07/22

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

| | | |
|---------------------------------------|------|---------------------------|
| Appearance: | OK | pale straw colored liquid |
| pH: (1% solution in DI water) | 12.1 | 11.5 - 13.0 |
| Surface Tension:(1% solution Dyne/cm) | 31.8 | 31.5 - 36.5 |
| Specific Gravity: | 1.27 | 1.26 - 1.34 |

Laura Falla, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.