



CERTIFICATE OF ANALYSIS
CITRANOX
FOR BECKMAN COULTER (BCI)

Date of Manufacture: 06/30/2020 to meet BC requirements
Manufactured by: Alconox, Inc.
Manufacturer Lot number: C9F0
BCI Lot number: C9F0
Manuf. Cat. No.: 1801B
BCI Part No: 81912
Shelf Life Expiration Date: 06/30/2022 to meet BC requirements
Storage Conditions: Closed upright and in a cool dry place (15 - 30 deg C or 60 - 90 deg F).

2022-06-30

REF 81912
LOT C9F0

+H62819121A

+S2063022C9F0A8



(01) 15099590223571
(11) 200630
(17) 220630
(10) C9F0

2020-06-30

Value	Result	Acceptance Range
Appearance:	OK	Clear pale olive colored liquid
pH: (1% solution in DI water)	2.8	2.0-3.0
Surface Tension: (1% solution Dyne/cm)	32.3	29.0 - 33.0
Specific Gravity	1.111	1.090 - 1.140
Number of Bar Code Labels Used for Labels	720	Not Applicable
Number of Bar Code Labels Used for Records	1	Not Applicable
Number of Bar Code Labels Destroyed	127	Not Applicable

Laura Falla

Laura Falla, QA