

CERTIFICATE OF ANALYSIS
CITRANOX
FOR BECKMAN COULTER (BCI)

Date of Manufacture: 09/30/2020 to meet BC requirements

Manufactured by: Alconox, Inc.

Manufacturer Lot number: C1G0

BCI Lot number: C9J0

Manuf. Cat. No.: 1801B

BCI Part No: 81912

Shelf Life Expiration Date: 09/30/2022 to meet BC requirements


Storage Conditions: Closed upright and in a cool dry place (15 - 30 deg C or 60 - 90 deg F).

<u>Value</u>	<u>-</u>	<u>Result</u>	<u>Acceptance Range</u>
Appearance:		OK	Clear pale olive colored liquid
pH: (1% solution in DI water)		2.9	2.0-3.0
Surface Tension: (1% solution Dyne/cm)		32.5	29.0 - 33.0
Specific Gravity		1.111	1.090 - 1.140
Number of Bar Code Labels Used for Labels		720	Not Applicable
Number of Bar Code Labels Used for Records		1	Not Applicable
Number of Bar Code Labels Destroyed		141	Not Applicable



Laura Falla, QA

REF 81912

 2022-09-30

LOT C9J0




+H628819121A



+\$\$2093022C9J0AF



(01)15099590223571
(11)200930
(17)220930
(10)C9J0

 2020-09-30

101550A

The Leader in Critical Cleaning
30 Glenn St. White Plains NY 10603 USA
Tel 914.948.4040 • Fax 914-948-4088
www.alconox.com • cleaning@alconox.com