



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 5/13/2021

## **CERTIFICATE OF ANALYSIS** **DETOJET**

**Manufactured by: Alconox Inc.**

Sub lot No: D1E1 D1E1F, D1E1G, D1E1H, D1E1J, D1E1K, D1E1L, D1E1M

Sub lot numbers of Master lot No: D1B0

Manuf. Cat. no.: 1601, 1601-1, 1605, 1615, 1630, 1632, 1632-1, 1655

Date Of Manufacture: 5/10/2021

Expires end of: 05/23

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

**ANALYSIS**

**RESULT**

**ACCEPTED RANGE**

Appearance:	OK	Clear, water white liquid
pH: (1% solution in DI water)	12.2	11.5 - 13.0
Foam (1% in DI water)	NIL	NIL
Specific Gravity:	1.34	1.27 - 1.34

Laura Falla, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*