



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 10/13/2021

## **CERTIFICATE OF ANALYSIS** **TERGAZYME**

**Manufactured by: Alconox Inc.**

Sub lot No: G1K1 G1K1A, G1K1B, G1K1C, G1K1D

Sub lot numbers of Master lot No: G1E0

Manuf. Cat. no.: 1301, 1303, 1304, 1304-1, 1325, 1350

Date Of Manufacture: 10/12/2021

Expires end of: 10/23

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

### **ANALYSIS**

### **RESULT**

### **ACCEPTED RANGE**

Appearance:	OK	White powder with cream colored flakes and small brown granules
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH:(1% solution in DI water)	9.2	9.0 - 10.0
Surface Tension:(1% solution Dyne/cm)	31.6	29.0 - 33.0

Laura Falla, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*