



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 12/15/2022

## CERTIFICATE OF ANALYSIS LIQUINOX

**Manufactured by: Alconox Inc.**

Sub lot No: L3M2 L3M2F, L3M2H, L3M2J, L3M2K, L3M2L, L3M2M

Sub lot numbers of Master lot No: L1B2

Manuf. Cat. no.: 1201, 1201-1, 1205, 1215, 1230, 1232, 1232-1, 1255

Date Of Manufacture: 12/11/2022

Expires end of: 12/24

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<u>ANALYSIS</u>	<u>RESULT</u>	<u>ACCEPTED RANGE</u>
Appearance:	OK	Clear golden liquid
pH: (as is)	8.1	8.0 - 9.5
Surface Tension:(1% solution Dyne/cm)	28.6	27.0 - 31.0
Specific Gravity:	1.08	1.06 - 1.09

Laura Falla, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*