



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 2/1/2023

CERTIFICATE OF ANALYSIS LIQUINOX

Manufactured by: Alconox Inc.

Sub lot No: L2A3 L2A3F, L2A3H, L2A3J, L2A3K, L2A3L, L2A3M

Sub lot numbers of Master lot No: L1B2

Manuf. Cat. no.: 1201, 1201-1, 1205, 1215, 1230, 1232, 1232-1, 1255

Date Of Manufacture: 1/30/2023

Expires end of: 01/25

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<u>ANALYSIS</u>	<u>RESULT</u>	<u>ACCEPTED RANGE</u>
Appearance:	OK	Clear golden liquid
pH: (as is)	8.3	8.0 - 9.5
Surface Tension:(1% solution Dyne/cm)	29.0	27.0 - 31.0
Specific Gravity:	1.07	1.06 - 1.09

Laura Falla, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.