



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 10/6/2023

**CERTIFICATE OF ANALYSIS**  
**SOLUJET**

**Manufactured by: Alconox Inc.**

Sub lot No: 10-03-23-02

Sub lot numbers of Master lot No: S1A2

Manuf. Cat. no.: 2101, 2101-1, 2105, 2115, 2130, 2155

Date Of Manufacture: 10/3/2023

Expires end of: 10/25

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<b><u>ANALYSIS</u></b>	<b><u>RESULT</u></b>	<b><u>ACCEPTED RANGE</u></b>
Appearance:	OK	pale straw colored liquid
pH: (1% solution in DI water)	11.9	11.5 - 13.0
Surface Tension:(1% solution Dyne/cm)	31.5	31.5 - 36.5
Specific Gravity:	1.27	1.26 - 1.34

Laura Falla, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*