



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 12/20/2023

CERTIFICATE OF ANALYSIS
CITRANOX

Manufactured by: Alconox Inc.

Sub lot No: 12-05-23-01

Sub lot numbers of Master lot No: C1A2

Manuf. Cat. no.: 1801, 1801-1, 1805, 1815, 1832, 1832-1, 1855

Date Of Manufacture: 12/5/2023

Expires end of: 12/25

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

| <u>ANALYSIS</u> | <u>RESULT</u> | <u>ACCEPTED RANGE</u> |
|---------------------------------------|----------------------|------------------------------|
| Appearance: | OK | Clear, olive colored liquid |
| pH: (1% solution in DI water) | 2.7 | 2.0 - 3.0 |
| Surface Tension:(1% solution Dyne/cm) | 31.2 | 29.0 - 33.0 |
| Specific Gravity: | 1.11 | 1.09 - 1.14 |

Laura Falla, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.