



TERGAZYME®

Enzyme-Active Powdered Detergent

Sample Trace Analysis Undiluted Concentrate

ELEMENT	DETECTED	DETECTION LIMIT
Chlorine	929 PPM	Less than 100 PPM
Fluorine	NONE	Less than 10 PPM
Bromine	35 PPM	Less than 10 PPM
Iodine	NONE	Less than 10 PPM
Boron	NONE	Less than 10 PPM
Lead	NONE	Less than 10 PPM
Zinc	NONE	Less than 10 PPM
Cadmium	NONE	Less than 50 PPM
Arsenic	NONE	Less than 50 PPM
Antimony	NONE	Less than 50 PPM
Bismuth	NONE	Less than 10 PPM
Tin	NONE	Less than 10 PPM
Sulfur	1210 PPM	Less than 100 PPM
Mercury	NONE	Less than 0.5 PPM
Sodium	An ingredient in TERGAZYME present at 20–30% by weight	
Phosphorus	An ingredient in TERGAZYME present at 7.6% by weight	

Note: TERGAZYME detergent is recommended for use at 1% dilutions.