



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 3/26/2024

CERTIFICATE OF ANALYSIS
ALCONOX

Manufactured by: Alconox Inc.

Sub lot No: 03-18-24-12

Sub lot numbers of Master lot No: A4A2

Manuf. Cat. no.: 1101, 1103, 1104, 1104-1, 1112, 1112-1, 1125, 1150

Date Of Manufacture: 3/18/2024

Expires end of: 03/26

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

| | | |
|---|------|--|
| Appearance: | OK | White powder with cream colored flakes |
| Solubility: (1% solution at 25 deg C, with agitation) | COMP | Complete |
| pH:(1% solution in DI water) | 9.3 | 9.0 - 10.0 |
| Surface Tension:(1% solution Dyne/cm) | 30.0 | 29.0 - 33.0 |

Laura Falla, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.