



Biodegradable. Concentrated. Interfering residue free. Alconox Inc. detergents set the standard for critical cleaning. Scientists, engineers and professionals specify Alconox Inc. detergents in their cleaning procedures for optimum economy, environmental safety and performance.

 = spray CIP, machine

 = manual cleaning


 = siphon pipet washer

 = soak, ultrasonic




USDA/FDA compliant. Dilute 1:100. pH 9.5

Size	Mfr. Cat No
4 lb box (1.8 kg)	1104-1
CS 9 x 4 lb (9 x 1.8 kg)	1104
25 lb box (11 kg)	1125
50 lb box (22.5 kg)	1150
100 lb drum (45 kg)	1101
300 lb drum (136 kg)	1103
50 x .5 oz packs (50 x 15g)	1112-1
CS 12 x 50 packs	1112



1 qt bottle (.95 L)	1232-1
CS 12 x 1 qt (12 x .95 L)	1232
1 gal bottle (3.8 L)	1201-1
CS 4 x 1 gal (4 x 3.8 L)	1201
5 gal jerrycan (19 L)	1205
15 gal drum (57 L)	1215
55 gal drum (209 L)	1255




and plastic. Surfactant rinse aids for no interfering residues. Has inhibitory residue test. USDA and FDA compliant. Dilute 1:100. pH 2.5

Size	Mfr. Cat No
1 gal bottle (3.8 L)	1801-1
CS 4 x 1 gal (4 x 3.8 L)	1801
5 gal jerrycan (19 L)	1805
15 gal drum (57 L)	1815
55 gal drum (209 L)	1855




and plastic. Surfactant rinse aids for no interfering residues. Has inhibitory residue test. USDA and FDA compliant. Dilute 1:100. pH 2.5

1 gal bottle	2301-1
CS 4 x 1 gal (4 x 3.8 L)	2301
5 gal jerrycan (19 L)	2305
55 gal drum (208 L)	2355



for noncorrosive, light-duty cleaning. No chelation or alkaline waste treat problems. Free rinsing, eliminating interfering residues. Dilute 3:100. pH 7.0


Size	Mfr. Cat No
1 gal bottle (3.8 L)	1901-1
CS 4 x 1 gal (4 x 3.8 L)	1901
15 gal drum (57 L)	1915
55 gal drum (209 L)	1955



Enzymatic cleaner for manual, soak, and ultrasonic cleaning of proteins, tissue, and body fluids from glass, metals, and plastic, with no interfering residues. Cleans filter membranes, fermenters and bioreactors. Has inhibitory residue test. FDA/USDA compliant. Dilute 1:100. pH 9.5

4 lb box (1.8 kg)	1304-1
CS 9 x 4 lb (9 x 1.8 kg)	1304
25 lb box (11 kg)	1325
50 lb box (22.5 kg)	1350
100 lb drum (45 kg)	1301
300 lb drum (136 kg)	1303

Available worldwide from leading laboratories, scientific, healthcare and industrial suppliers.



CITRAJET® Low-Foam Liquid Acid Cleaner/Rinse



Liquid acid cleaner for machine, spray-CIP, manual, and ultrasonic cleaning or acid neutralizing. Corrosion-inhibited for glassware, metals, and plastic. Surfactant rinse for no interfering residues. Has inhibitory residue test. FDA compliant. Dilute 1:100. pH 2.5

Size	Mfr. Cat No
1 gal bottle (3.8 L)	2001-1
CS 4 x 1 gal (4 x 3.8 L)	2001
5 gal jerrycan (19 L)	2005
15 gal drum (57 L)	2015
55 gal drum (209 L)	2055



DETOJET® Low-Foaming Liquid Detergent



Heavy-duty alkaline degreaser. Nonionic liquid for labware and machine washers, manual, and ultrasonic cleaning. Corrosion-inhibited for glass, metals, and plastic. No interfering residues. Has inhibitory residue test. USDA/FDA compliant. Dilute 1:200. pH 12

1 gal bottle (3.8 L)	1601-1
CS 4 x 1 gal (4 x 3.8 L)	1601
5 gal can (19 L) - Air OK	1605
15 gal drum (57 L)	1615
30 gal drum (114 L)	1630
55 gal drum (209 L)	1655



ALCOJET® Low-Foaming Powdered Detergent



General-purpose nonionic powder for labware and machine washers, manual, and ultrasonic cleaning. Corrosion-inhibited for glassware, metals, and plastic. No interfering residues. Has inhibitory residue test. USDA authorized. Dilute 1:200. pH 11.5

Size	Mfr. Cat No
4 lb box (1.8 kg)	1404-1
CS 9 x 4 lb (9 x 1.8 kg)	1404
25 lb box (11 kg)	1425
50 lb box (22.5 kg)	1450
100 lb drum (45 kg)	1401
300 lb drum (136 kg)	1403



ALNOCHROMIX™ Oxidizing Acid Additive for Glass Cleaning



Concentrated inorganic persulfate based additive to sulfuric acid to make a metal-free and surfactant-free glassware cleaner to be used in tissue culture, trace metal, and enzyme kinetics research, etc. Nonmetallic, non-carcinogenic, and easier to dispose replacement for chromic acid.

1 box 10 x 3.1 oz (87 g) pouches	2510
4 boxes of 10 x 3.1 oz (87 g) pouches	2540
1 pair silver gloves	2521
3 safety vent caps	2513
Kit 3 vent cap+3 pair gloves	2533



TERGAJET® Low-Foaming Phosphate-Free Powder



Phosphate-free, nonionic powder for machine washers, hand, and ultrasonic cleaning. Corrosion-inhibited for labware. For water, microbiology, and environmental labs. No interfering residues. Has inhibitory residue test. Dilute 1:200. pH 11.5

Size	Mfr. Cat No
4 lb box (1.8 kg)	2204-1
CS 9 x 4 lb (9 x 1.8 kg)	2204
25 lb box (11 kg)	2225
50 lb box (22.5 kg)	2250
100 lb drum (45 kg)	2201
300 lb drum (136 kg)	2203



SOLUJET® Low-Foaming Phosphate-Free Liquid



Nonionic liquid for machine, spray-CIP, manual, soak, and ultrasonic cleaning. High emulsifying and penetrating formula. Corrosion-inhibited for glass, metals, and plastic. No interfering residues. Has inhibitory residue test. FDA compliant. Dilute 1:200. pH 12

1 gal bottle (3.8 L)	2101-1
CS 4 x 1 gal (4 x 3.8 L)	2101
5 gal can (19 L) - Air OK	2105
15 gal drum (57 L)	2115
55 gal drum (209 L)	2155



DETERGENT 8® Low-Foaming Ion-Free Liquid Detergent

Nonionic detergent for PC board washing machines, manual, and ultrasonic cleaning. Ion- and chelate-free for no conductive residues on electronics and no chelate interference in waste. Phosphate-free. Dilute 3:100. pH 11

Size	Mfr. Cat No
1 gal bottle (3.8 L)	1701-1
CS 4 x 1 gal (4 x 3.8 L)	1701
5 gal jerrycan (19 L) - Air OK	1705
15 gal drum (57 L)	1715
55 gal drum (209 L)	1755



KEYLAJET® Low-Foaming High Alkaline Liquid



Corrosion inhibited formulation recommended for hard metals, steels, stainless steel, porcelain, ceramic, rigid plastics, cement, and elastomers. Free rinsing to give you reliable results and no interfering residues. Dilute 1:200. pH 13+

5 gal can (19 L) - Air OK	2405
55 gal drum (208 L)	2455
265 gal tote (1003 L)	24265



ALCOTABS® Tablet Pipet Detergent



Effervescent detergent tablet for maximum cleaning in siphon pipet washers. Slow release bubbles and free rinsing give no interfering residues. Has inhibitory residue test. Use 2 tabs/wash. pH 6.5

100 tablet bottle	1500-1
CS 6 x 100 bottles	1500